



Partisep™ – ES Industries has developed a series of chemistries to closely match the performance of GE Whatman Partisil® columns. Partisep materials are available in both the standard 10 micron particle size as well as a 5 micron particle size for both shorter run times and higher efficiencies. Most methodologies on Partisil HPLC columns can be transferred to Partisep products without modification, including USP applications. The following Partisep phases are available:

- Partisep ODS
- Partisep ODS2
- Partisep ODS3
- Partisep C8
- Partisep Silica
- Partisep PAC
- Partisep SAX
- Partisep SCX

Partisep™ ODS 5u 85A				
	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-ODS	114271-PSP-ODS	11d271-PSP-ODS	112271-PSP-ODS
10cm	125271-PSP-ODS	124271-PSP-ODS	12d271-PSP-ODS	122271-PSP-ODS
15cm	135271-PSP-ODS	134271-PSP-ODS	13d271-PSP-ODS	132271-PSP-ODS
25cm	155271-PSP-ODS	154271-PSP-ODS	15d271-PSP-ODS	152271-PSP-ODS

Partisep™ ODS2 5u 85A				
	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-ODS2	114271-PSP-ODS2	11d271-PSP-ODS2	112271-PSP-ODS2
10cm	125271-PSP-ODS2	124271-PSP-ODS2	12d271-PSP-ODS2	122271-PSP-ODS2
15cm	135271-PSP-ODS2	134271-PSP-ODS2	13d271-PSP-ODS2	132271-PSP-ODS2
25cm	155271-PSP-ODS2	154271-PSP-ODS2	15d271-PSP-ODS2	152271-PSP-ODS2

Partisep™ ODS3 5u 85A				
	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-ODS3	114271-PSP-ODS3	11d271-PSP-ODS3	112271-PSP-ODS3
10cm	125271-PSP-ODS3	124271-PSP-ODS3	12d271-PSP-ODS3	122271-PSP-ODS3
15cm	135271-PSP-ODS3	134271-PSP-ODS3	13d271-PSP-ODS3	132271-PSP-ODS3
25cm	155271-PSP-ODS3	154271-PSP-ODS3	15d271-PSP-ODS3	152271-PSP-ODS3

Partisep™ C8 5u 85A

	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-C8	114271-PSP-C8	11d271-PSP-C8	112271-PSP-C8
10cm	125271-PSP-C8	124271-PSP-C8	12d271-PSP-C8	122271-PSP-C8
15cm	135271-PSP-C8	134271-PSP-C8	13d271-PSP-C8	132271-PSP-C8
25cm	155271-PSP-C8	154271-PSP-C8	15d271-PSP-C8	152271-PSP-C8

Partisep™ Silica 5u 85A

	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-SI	114271-PSP-SI	11d271-PSP-SI	112271-PSP-SI
10cm	125271-PSP-SI	124271-PSP-SI	12d271-PSP-SI	122271-PSP-SI
15cm	135271-PSP-SI	134271-PSP-SI	13d271-PSP-SI	132271-PSP-SI
25cm	155271-PSP-SI	154271-PSP-SI	15d271-PSP-SI	152271-PSP-SI

Partisep™ PAC 5u 85A

	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-PAC	114271-PSP-PAC	11d271-PSP-PAC	112271-PSP-PAC
10cm	125271-PSP-PAC	124271-PSP-PAC	12d271-PSP-PAC	122271-PSP-PAC
15cm	135271-PSP-PAC	134271-PSP-PAC	13d271-PSP-PAC	132271-PSP-PAC
25cm	155271-PSP-PAC	154271-PSP-PAC	15d271-PSP-PAC	152271-PSP-PAC

Partisep™ SAX 5u 85A

	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-SAX	114271-PSP-SAX	11d271-PSP-SAX	112271-PSP-SAX
10cm	125271-PSP-SAX	124271-PSP-SAX	12d271-PSP-SAX	122271-PSP-SAX
15cm	135271-PSP-SAX	134271-PSP-SAX	13d271-PSP-SAX	132271-PSP-SAX
25cm	155271-PSP-SAX	154271-PSP-SAX	15d271-PSP-SAX	152271-PSP-SAX

Partisep™ SCX 5u 85A

	<u>4.6mm (ID)</u>	<u>4.0mm (ID)</u>	<u>3.2mm (ID)</u>	<u>2.1mm (ID)</u>
5cm	115271-PSP-SCX	114271-PSP-SCX	11d271-PSP-SCX	112271-PSP-SCX
10cm	125271-PSP-SCX	124271-PSP-SCX	12d271-PSP-SCX	122271-PSP-SCX
15cm	135271-PSP-SCX	134271-PSP-SCX	13d271-PSP-SCX	132271-PSP-SCX
25cm	155271-PSP-SCX	154271-PSP-SCX	15d271-PSP-SCX	152271-PSP-SCX

ES INDUSTRIES
701 South Route 73
West Berlin NJ 08091

Toll-Free: 1-800-356-6140
Telephone: 856-753-8400
Fax: 856-753-8484

email: esindustries@msn.com website: www.esind.com